INSTRUCTIONS FOR USE
DynamX Novolimus Eluting Coronary Bioadaptor System

TABLE OF CONTENTS

ENGLISH (EN).....	2
DEUTSCH (DE).....	4
DANSK (DA)	7
ESPAÑOL (ES).....	10
SVENSKA (SV).....	12
简体中文 (ZH-S).....	15
繁體中文 (ZH-T).....	17
العربية (AR).....	19
FRANÇAIS (FR)	21
ITALIANO (IT)	24
PORTUGUÊS (PORTUGAL) (PT).....	26
ΕΛΛΗΝΙΚΑ (EL)	29
NEDERLANDS (NL).....	32
SUOMI (FI).....	34
NORSK (NO).....	37
POLSKI (PL).....	39
SLOVENČINA (SK)	42
ČEŠTINA (CS)	44
TÜRKÇE (TR).....	47
Русский (RU)	49
עברית (HE).....	52
PORTUGUÊS (BRASIL) (PT-BR)	54
ROMÂNĂ (RO).....	56
HRVATSKI (HR)	59
EESTI (ET).....	61
LATVISKI (LV)	64
LİETUVIŲ K. (LT)	66
БЪЛГАРСКИ (BG)	68

Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.
 For export only - not for sale in the USA

This device should only be used by physicians trained in angiography and percutaneous transluminal coronary angioplasty (PTCA).

1.0 DEVICE DESCRIPTION

The DynamX Novolimus Eluting Coronary Bioadaptor System (DynamX NECBS) is comprised of two key components: a metallic bioadaptor coated with a biodegradable polymer and drug (Novolimus) and the delivery system.

The DynamX NECBS are comprised of a pre-mounted balloon expandable cobalt chromium alloy bioadaptor. The bioadaptor pattern includes expansion segments evenly spaced along the length and is coated with a proprietary bioresorbable polylactide based polymer. Approximately 5 mcg of Novolimus per mm of bioadaptor length is delivered via the biodegradable polymer top coat which allows sustained release of the majority of the drug in approximately 4 weeks.

The DynamX delivery system incorporates a rapid-exchange delivery system that includes a nylon-blend balloon and distal shaft and a stainless steel hypotube in the proximal shaft. The total working length of the delivery system is 140 cm. There are two radiopaque markers located underneath the balloon to aid in visualization and to facilitate accurate bioadaptor placement. Additionally, there are two proximal delivery system shaft markers (90 cm and 100 cm from the distal tip) that indicate the relative position of the delivery system to the end of a brachial or femoral guiding catheter tip. The delivery system is compatible with 0.014 inch (0.36 mm) guide wires and 5 French guiding catheters (minimum inner diameter of 0.058 inches, or 1.5 mm).

The DynamX NECBS label specifications are included in summary format in Table 1.0 below.

Table 1.0 DynamX NECBS Label Specifications

Bioadaptor Length(s) (Nominal)	Drug Dose (Nominal)	Minimum Guiding Catheter Compatibility (ID)	In Vitro Bioadaptor Rated Burst Pressure (RBP)*
14 mm	65 mcg	5F (0.058")	16 atm
18 mm	85 mcg		
23 mm	105 mcg		
28 mm	125 mcg		
32 mm	145 mcg		
38 mm	170 mcg		

*See package label for product compliance information.

2.0 HOW SUPPLIED

- Sterile:** This device is sterilized with electron beam radiation and is non-pyrogenic. **Do not use the device if the package is opened or damaged.**
- Contents:** One DynamX Novolimus Eluting Coronary Bioadaptor System
 One flushing needle
 One Instructions for Use document
- Storage:** Store at ≤ 25°C.

3.0 INTENDED USE AND INDICATIONS FOR USE

3.1 Intended Use

The DynamX Bioadaptor is a permanent implant incorporating the ancillary drug substance, Novolimus, and is intended to improve coronary luminal diameter in patients with symptomatic ischemic heart disease.

3.2 Indications for Use

The DynamX NECBS is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease due to discrete de novo native coronary artery lesions with reference vessel diameters between 2.25 - 3.5 mm and ≤ 34 mm in length. The bioadaptor is intended as a permanent implant.

4.0 CONTRAINDICATIONS

The DynamX NECBS is contraindicated for use in:

- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon
- Persons allergic to cobalt chromium alloy (including the major elements cobalt, chromium, tungsten, nickel), polylactide, polyglycolide, or polycaprolactone based polymers, or drugs such as Novolimus, Rapamycin or structurally related compounds may suffer allergic reaction to this implant

5.0 WARNINGS AND PRECAUTIONS

WARNINGS

- Judicious selection of patients is necessary since the use of this device carries the associated risk of thrombosis, vascular complications and/or bleeding events.
- Persons who are unable to tolerate long-term dual antiplatelet therapy for a minimum of 12 months should not be considered for implantation of the DynamX NECBS.
- Only physicians who have received appropriate training should perform implantation of the bioadaptor.
- Bioadaptor placement should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- Subsequent bioadaptor blockage may require repeat dilatation of the arterial segment containing the bioadaptor. The long term outcome following repeat dilatation of endothelialized bioadaptors is not well characterized.

- To avoid the possibility of dissimilar metal corrosion, do not implant devices of different materials in tandem where overlap or contact is possible.
- Consideration should be given to the risks and benefits of use in patients with history of severe reaction to contrast agents.
- The safety and effectiveness of using mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters in conjunction with DynamX NECBS implantation have not been established.

PRECAUTIONS

Intended for single use only. Reuse, reprocess or resterilization of the device may result in loss of drug coating efficacy and/or loss of device integrity and/or contamination which may result in repeat intervention and/or patient injury, illness, or death.

5.1 Bioadaptor Handling

- For single use only.** Do not resterilize or reuse this device. Note product “Use by” date on the label.
- Do not remove the Bioadaptor from the Delivery System** as removal may damage the bioadaptor and/or lead to bioadaptor embolization. The DynamX NECBS is intended to perform as a system. The bioadaptor is not designed to be reloaded onto another delivery device. Attempts to reload the bioadaptor may damage the bioadaptor and/or lead to embolization.
- Special care must be taken not to handle or in any way disrupt the bioadaptor position on the Delivery System. This is most important during catheter removal from the packaging, placement over the guide wire and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- Do not manipulate (e.g. roll) the bioadaptor with your fingers, as this action may loosen the bioadaptor from the delivery system.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon, as this may cause uneven expansion and difficulty in deployment of the bioadaptor.

5.2 Bioadaptor Placement

- Predilatation of the lesion is required for safe bioadaptor delivery.
- Do not prepare or pre-inflate the bioadaptor delivery system prior to bioadaptor deployment** other than as directed. Use balloon purging technique described in Preparation of Delivery System section.
- Implanting a bioadaptor may lead to dissection of the vessel proximal or distal to the bioadaptor-implanted portion and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilatation or placement of additional bioadaptors). If additional bioadaptor placement is required, bioadaptor materials of similar composition should be used.
- When treating multiple lesions, place the bioadaptor in the distal lesion prior to placement of a bioadaptor in the proximal lesion to reduce the chance of dislodging the proximal bioadaptor.
- Do not expand the bioadaptor if it is not properly positioned in the vessel (See Bioadaptor/Delivery System Removal Precautions).
- Bioadaptor placement may compromise side branch patency.
- Do not exceed the Rated Burst Pressure as indicated on the product label.** Balloon pressures should be monitored during inflation. Use of pressures higher than those specified on the product label may result in balloon rupture and/or potential vessel damage.
- Should any resistance be felt at any time during either lesion access or removal of the delivery system pre-bioadaptor implantation, the system should be removed as a single unit per the Bioadaptor/Delivery System Removal Precautions.
- Bioadaptor retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site.
- Additional expansion of a deployed bioadaptor may cause a flow limiting dissection. This may be treated by implantation of another bioadaptor. When multiple bioadaptors are used, the ends should overlap slightly.
- The vascular effect of overlapping bioadaptors has not been clinically tested with this device and should be employed only as medically necessary.
- Ensure the bioadaptor is not under dilated. If the deployed bioadaptor is not fully apposed to the vessel wall, the bioadaptor may be expanded further with a high pressure non-compliant balloon that is slightly shorter than the bioadaptor. The edges of the balloon should not extend outside of the bioadaptor-implanted region.

5.3 Bioadaptor / Delivery System Removal

Should **any resistance** be felt **at any time** during either lesion access or removal of the Delivery System before bioadaptor implantation, the Delivery System and guiding catheter should be removed as a single unit.

When removing the delivery system as a single unit:

- Do not retract (withdraw) the Delivery System into the guiding catheter.
- Position the proximal balloon marker just distal to the tip of the guiding catheter.
- Advance the guidewire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter; then remove the guiding catheter and delivery system as a **single unit**.

Failure to follow these steps and/or applying excessive force to the Delivery System can potentially result in loss or damage to the bioadaptor and/or delivery system components.

If it is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

5.4 Post-Implant Procedure

When crossing a newly deployed bioadaptor with a guide wire, balloon or delivery system, exercise care to avoid disrupting the bioadaptor geometry.

5.4.1 Magnetic Resonance Imaging (MRI)

The DynamX Novolimus Eluting Coronary Bioadaptor has not been tested for safety in the MRI environment. Therefore, MRI scans should not be performed on patients post-implantation until the bioadaptor has completely endothelialized to minimize the potential for migration. For a conventional uncoated CoCr cobalt chromium stent, this period is usually considered to be 8 weeks. This time may be delayed due to the presence of the drug and polymer coating. This device has not been evaluated for heating in the MRI environment. The effect of heating in the MRI environment on the drug and polymer coating is not known. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the bioadaptor.

6.0 DRUG INTERACTION

While no specific clinical data are available, drugs like tacrolimus which act through the same binding protein (FKBP) may interfere with the efficacy of Novolimus or Rapamycin-type drugs. Drug interaction studies have not been performed. Rapamycin-type drugs are metabolized by CYP3A4. Strong inhibitors of CYP3A4 (e.g. ketoconazol) might cause increased Novolimus exposures to levels associated with systemic effects, especially if multiple bioadaptors are deployed. Systemic exposure of Novolimus or Rapamycin-type drugs should also be taken into consideration if the patient is treated concomitantly with systemic immunosuppressive therapy. Grapefruit may potentially interfere with the metabolism of Novolimus.

7.0 POTENTIAL ADVERSE EVENTS

Adverse events which may be associated with the implantation of a coronary bioadaptor or PTCA procedure include (in alphabetical order):

- Abrupt Closure
- Allergic reaction (hypersensitivity to polylactide, polyglycolide, or polycaprolactone based polymers, or cobalt chromium alloy)
- Aneurysm, pseudoaneurysm
- Angina
- Arrhythmias
- Arteriovenous fistula
- Bioadaptor migration, embolization or misplacement
- Cardiac tamponade
- Damage to the vessel requiring CABG
- Death
- Dissection, perforation or rupture of the coronary artery
- Drug reactions to antiplatelet/anticoagulation agents or contrast media
- Embolism (air, tissue, device or thrombus)
- Failure to deliver the device to the intended site
- Hemorrhage
- Hypotension/hypertension
- Infection and pain at the insertion site
- Myocardial ischemia and/or infarction
- Occlusion of the artery
- Peripheral ischemia
- Restenosis of the bioadapted artery
- Stroke/transient ischemic attack
- Thrombosis (acute, subacute or late)
- Ventricular fibrillation
- Vessel spasm

The following additional side effects/complications may be associated with, but not limited to the use of Novolimus or Rapamycin-type drugs:

- Acne
- Diarrhea or constipation
- Headache
- Increased blood pressure
- Increased cholesterol or triglyceride levels
- Insomnia
- Nausea
- Rash
- Sore or weak muscles or joints
- Tremor
- Upper respiratory or urinary tract infections
- Water retention

8.0 PATIENT SELECTION AND TREATMENT

8.1 Individualization of Treatment

The risks and benefits should be carefully considered for each patient before use of the DynamX NECBS. Patient selection factors to be assessed should include a judgment regarding risk of anticoagulation therapy. Persons who are unable to tolerate long-term dual antiplatelet therapy for a minimum of 12 months should not be considered for implantation of the DynamX NECBS. Bioadaptor implantation should be generally avoided in those patients at heightened risk of bleeding (e.g., those patients with recently active gastritis or peptic ulcer disease) in which anticoagulation therapy would be contraindicated.

Co-morbidities that increase the risk of poor initial results or the risk of emergency bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed.

Thrombosis following bioadaptor implantation is affected by several baseline angiographic and procedural factors. These include vessel diameters less than 2.0 mm, intra-procedural thrombosis, poor distal flow, and/or dissection following bioadaptor implantation. In patients that have undergone coronary bioadaptor implantation, the persistence of thrombus or dissection is considered a marker for subsequent thrombotic occlusion. These patients should be monitored very carefully during the first month after bioadaptor implantation.

Use outside of the specified Indications for Use is prohibited. The use of drug-eluting stents/devices in patients and lesions outside the labeled indications may have an increased risk of adverse events, including device thrombosis, device embolization, myocardial infarction, or death.

8.2 Use in Specific Patient Populations

The safety and effectiveness of the DynamX NECBS have not been established in patients with:

- Unresolved thrombus at the lesion site.
- A coronary artery reference vessel diameter < 2.25 mm.
- A lesion located in the left main coronary artery, ostial lesions or lesions located at a bifurcation.
- Diffuse disease, or poor outflow distal to the identified lesions.
- > 2 overlapping bioadaptors (due to risk of thrombosis and restenosis).

The safety and effectiveness of using mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters to treat in-bioadaptor restenosis have not been established.

8.3 Pregnancy

Pregnancy Category C: There are no adequate and well controlled studies of Novolimus or DynamX bioadaptors in pregnant women.

Effective contraception should be initiated before implanting a DynamX bioadaptor and for 12 weeks after implantation. The DynamX bioadaptor should be used during pregnancy only if the potential benefit outweighs the risk to the embryo or fetus.

8.4 Lactation

It is not known whether Novolimus is excreted in human milk. The pharmacokinetic and safety profiles of Novolimus in infants are not known. Because similar drugs are excreted in human milk and because of the potential for adverse reactions in nursing infants from Novolimus, a decision should be made whether to discontinue nursing or to implant the bioadaptor, taking into account the importance of the bioadaptor to the mother.

9.0 CLINICIAN USE INFORMATION

9.1 Materials Required

- Appropriate guiding catheter(s) (5F or 0.058" ID minimum)
- 2-3 syringes
- Heparinized normal sterile saline
- 0.014 inch (0.36 mm)(maximum OD) x 175 cm (minimum length) guide wire
- Rotating hemostatic valve
- 60% contrast diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque device
- Guide wire introducer

9.2 Packaging Removal

9.2.1 Carefully inspect the packaging pouch for damage. **Do not use if the package has been damaged or opened.**

9.2.2 Peel open the pouch using aseptic technique to reveal the sterile catheter.

9.2.3 Pass or drop the sterile catheter into the sterile field using aseptic technique.

9.3 Inspection Prior to Use

9.3.1 Remove the bioadaptor system from the hoop.

9.3.2 Prior to use, carefully inspect the catheter for kinks, bends or other damage. Do not use if damage is noted or if the device is inadvertently contaminated.

9.3.3 Verify that the bioadaptor is located between the radiopaque markers on the balloon.

9.4 Preparation of Delivery System

9.4.1 Prepare the guiding catheter, guide wire and inflation device according to the manufacturer's instructions.

9.4.2 Carefully remove the protective sheath covering the bioadaptor/balloon segment.

9.4.3 Flush the bioadaptor system with heparinized saline using the flushing needle.

Caution: Avoid manipulation of the bioadaptor during flushing of the guidewire lumen, as this may disrupt the placement of the bioadaptor on the balloon.

9.4.4 Fill a 20 cc syringe with 5 cc of contrast/heparinized normal saline mixture (1:1).

9.4.5 Attach to the delivery system proximal hub and apply negative pressure for 10 seconds.

9.4.6 With the syringe tip pointing downwards, slowly release pressure to neutral.

9.4.7 Attach a prepared inflation device to the catheter hub, leave on neutral pressure.

Caution: Do not apply negative pressure on inflation device during delivery of the bioadaptor to the lesion site.

9.5 Delivery Procedure

- 9.5.1 Prepare vascular access site according to standard practice.
- 9.5.2 Pre-dilate the lesion site using a balloon that is shorter than the bioadaptor length.
- 9.5.3 While maintaining neutral pressure on the inflation device, backload the delivery system onto the guidewire.
- 9.5.4 Fully open the rotating hemostatic valve to allow easy passage of the bioadaptor and delivery system.

Note: If resistance is felt at any time, determine the cause of resistance before proceeding and as necessary follow the Bioadaptor Removal Precaution instructions.

- 9.5.5 Tighten hemostatic valve and advance the bioadaptor system over the guidewire to the lesion site while maintaining stable guide catheter and guidewire positioning.

9.6 Bioadaptor Deployment

- 9.6.1 Prior to bioadaptor expansion, reconfirm the location of the bioadaptor in relation to the balloon markers using high-resolution fluoroscopy.
- 9.6.2 Under fluoroscopic guidance, inflate the bioadaptor system to the nominal pressure or until the bioadaptor appears fully expanded. Maintain pressure for 15 – 30 seconds. Optimal expansion requires the bioadaptor to be in full contact with the vessel wall, with the inner diameter of the bioadaptor matching the reference vessel diameter.

Note: Do not exceed the Rated Burst Pressure. Do not expand the bioadaptor beyond 0.5 mm of the nominal diameter.

- 9.6.3 If further dilatation of the deployed bioadaptor is required, a larger diameter, non-compliant balloon of a length shorter than the deployed bioadaptor may be used. Avoid disruption of the bioadaptor when crossing the newly deployed bioadaptor with the balloon catheter.

9.7 Removal Procedure

- 9.7.1 Ensure the balloon is fully deflated before attempting to withdraw the delivery system. Depending on device size and contrast brand and dilution, this may take up to 30 seconds.
- 9.7.2 Withdraw the delivery system while maintaining guidewire placement.
- 9.7.3 Repeat angiography to assess the bioadaptor-implanted area. If adequate expansion has not been achieved follow the instructions under Bioadaptor Deployment. **Ensure the bioadaptor is not under dilated.**

9.8 Bioadaptor/System Removal – Precautions

Should **any resistance** be felt **at any time** during either lesion access or removal of the bioadaptor delivery system pre-implantation, **remove the delivery system as a single unit.**

- Position the proximal balloon marker just distal to the tip of the guiding catheter.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the bioadaptor delivery system to the guiding catheter, and then remove the guiding catheter and bioadaptor system as a single unit.

Failure to follow these steps and/or applying excessive force to the bioadaptor delivery system during removal may result in the loss or damage to the bioadaptor and/or delivery system components.

10.0 PACKAGING INFORMATION

A package contains one DynamX Bioadaptor pre-loaded onto the Delivery System and one flushing needle. An Instructions for Use document is included in the package.

The DynamX NECBS is supplied Sterile using electron beam radiation and non-pyrogenic in unopened, undamaged packages.

Intended for single use only. Do not reuse, reprocess or resterilize.

Store at ≤ 25°C. Use by the “Use by” Date noted on the package.

CAUTION: DO NOT USE IF THERE IS ANY DAMAGE TO THE PACKAGE

11.0 PATENTS

U.S. and foreign patents pending

12.0 DISCLAIMER OF WARRANTY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. ELIXIR MEDICAL SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND ELIXIR MEDICAL TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the Disclaimer of Warranty shall not be affected.

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Elixir Medical Corporation
Acorn Regulatory Consultancy Services Ltd.