English [EN] Instructions for use - TrapIT Trapping Balloon Catheter	1
Español [ES] Instrucciones de uso - Catéter TrapIT Trapping Balloon Catheter	3
Français [FR] Mode d'emploi - TrapIT Trapping Balloon Catheter	4
Deutsch [DE] Gebrauchsanweisung - TrapIT Trapping Balloon Catheter	6
Italiano [IT] Istruzioni per l'uso - TrapIT Trapping Balloon Catheter	7
Türkçe [TR] Kullanım talimatları - TrapIT Trapping Balloon Kateteri	9
Slovensky [SK] Návod na použitie - Zachytávací balónový katéter TrapIT	10
Português [PT] Instruções de utilização - Cateter balão TrapIT	12
Polski [PL] Instrukcja użytkowania - Cewnik balonowy TrapIT Trapping	13
yesky [CS] Návod k použití - Balónkový katétr TrapIT	15
Dansk [DA] Brugervejledning - TrapIT Trapping-balloonkateter	16
Suomi [FI] Käyttöohjeet - Pysäytävä TrapIT-pallokatetri	18
Norsk [NO] Bruksanvisning - TrapIT-innfangningsballonkateter	19
Svenska [SV] Bruksanvisning - TrapIT Läsballonkateter	20
[BG]	22
Eesti [ET] Kasutusjuhend - TrapIT lukustuv balloonkateter	24
Magyar [HU] Használati utasítás - TrapIT befogó ballonkatéter	25
Latviešu [LV] Lietošanas pam. c ba - TrapIT nosiešanas balonkatērs	26
Lietuvi k. [LT] Naudojimo instrukcijos - „TrapIT“ rezekcinis balioninis kateteris	28
Română[RO] Instruc iuni de utilizare - Cateter cu balon de fixare TrapIT	29
Slovenščina [SL] Navodila za uporabo - Prestrezní balonski kateter TrapIT	31

koččowki dystalnej (TRP10015) [CS] Značka hloubky 1015mm od distálního hrotu (TRP10015) [DA] Dybdemarkering 1015mm fra den distale spids (TRP10015) [FI] Syvyysmerkki 1015mm:n kohdalla distaalisesta kärjestä (TRP10015) [NO] Dybdemarkering ved 1015mm fra den distale spissen (TRP10015) [SV] Djupmarkering vid 1015mm från den distala spetsen (TRP10015) [BG] 1015 mm (TRP10015) [ET] Sügavusmärgistus 1015 mm kaugusel distaalsest otsakust (TRP10015) [HU] Mélységjelölés a distális hegytől 1015 mm-re (TRP10015) [LV] Dzi uma atz me 1015 mm att lum no dist l gala (TRP10015) [LT] Gylio žyma 1015 mm nuo distalinio galo (TRP10015) [RO] Marcaj de adăncime la 1015 mm de la vârful distal (TRP10015) [SL] Oznaka globine na 1015 mm od distalne konice (TRP10015)
7. [EN] Hard stop [ES] Tope [FR] Butée dure [DE] Anschlagpuffer [IT] Punto di arresto [TR] Sert durdurucu [SK] Tvrdý doraz [PT] Paragem rígida [PL] Ogranicznik [SV] Pevná zarážka [DA] Hårdt stop [FI] Pysäytyskohta [NO] Hard stopper [SV] Stoppanordnin [BG] [ET] Peataja [HU] Kemény ütközés [LV] Piln gais atdurs [LT] Galutinis stabilizikis [RO] Blocare [SL] Trda zaustavitel
8. [EN] Strain relief [ES] Aliviador de tensión [FR] Détendeur [DE] Zugentlastung [IT] Dispositivo antitensione [TR] Gerilim azaltma [SK] Kompenzátor napnutia [PT] Dispositivo de alívio de tensão [PL] Zaczisk oddziałający [CS] Prvek pro odlehbení tahu [DA] Aflastning [FI] Vahvike [NO] Strekkavlastning [SV] Brytskyd [BG] [ET] Pingevastabi [HU] Feszültségmentesítő [LV] Atslodzes uzmaiva [TR] tampos maşinimas [RO] Eliberarea presiunii [SL] Natezna razbremenitev
9. [EN] Hub [ES] Conector del catéter [FR] Embase [DE] Anschlussstück [IT] Connettore [TR] Göbek [SK] Náboj [PT] Conector [PL] Nasadka [CS] Rozbojováb [DA] Hub [FI] Kanta [NO] Nav [SV] Kopplingsstyck [BG] [ET] Konnektor [HU] Elcsőzt [LV] Rumba [LT] vor [RO] Racord [SL] Nastavek

English [EN] Instructions for use - TrapIT Trapping Balloon Catheter

CAUTION

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

Device Description

The TrapIT Trapping Balloon Catheters have an integrated shaft system and a balloon near the distal tip. The shaft has a lumen which is used for inflation of the balloon with contrast medium. As aid in positioning the balloon in a guide catheter a tactile feedback stop is integrated in the shaft. The shaft of the 100 cm compatible TrapIT has a shaft depth marking that aid in gauging balloon catheter position relative to the guiding catheter tip when used in 90 cm guide catheters. The design of this balloon catheter does not incorporate a lumen for distal dye injections and distal pressure measurements.

STERILE. Sterilized with Gamma Radiation. Non-pyrogenic. For one time use only. Do not re-sterilize.

Eudamed and SSCP

A Summary of the Safety and Clinical Performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the BasicUDI-DI: 871848174TRP0001XZ.

The TrapIT is available in the following size

Model	TRP10015	TRP9015
Balloon length	15mm	15mm
Balloon diameter at nominal	2.5mm	2.5mm
Effective length	1115mm	1015mm
Depth marker on shaft	1015mm	-
Radiopaque tip length	4mm	4mm
Nominal pressure	8 atm (811 kPa)	8 atm (811 kPa)
Rated burst pressure	12 atm (1216 kPa)	12 atm (1216 kPa)

Intended purpose

The TrapIT is a trapping balloon catheter intended to facilitate interventional device exchange while maintaining guide wire position in patients undergoing PCI procedures.

Indications for use

The TrapIT is indicated for patients suffering from coronary or peripheral artery disease.

Contraindications


None.

Warnings

- For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death.
- Use the catheter before the "Use By" date specified on the package.
- Do not use if the package is open or damaged.
- Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of reuse, cleaning, and / or re-sterilization may compromise the structural integrity of the device and / or materials, leading to contamination due

Device Name	Model	Compatible Guiding catheter length	UDI DI
TrapIT 90	TRP9015	90 cm	08718481740262
TrapIT 100	TRP10015	90 and 100 cm	08718481740279

Note: Not all device models are available in every country.

1. [EN] Distal tip [ES] Punta distal [FR] Extrémité distale [DE] Distale Spitze [IT] Punta distale [TR] Distal uç [SK] Distální koniec [PT] Ponta distal [PL] Kończówka dystalna [CS] Distální hrot [DA] Distal spids [FI] Distaalinen kärki [NO] Distal spiss [SV] Distal spet [BG] [ET] Distalne otsak [HU] Disztális hegy [LV] Dist lais gals [LT] Distalinis galas [RO] Vârful distal [SL] Distalna konica
2. [EN] Radiopaque marker [ES] Marcador radiopaco [FR] Marqueur radio-opaque [DE] Röntgendichte Markierung [IT] Indicatore radiopaco [TR] Radyopak iüaretleiyici [SK] RTG kontrastná značka [PT] Marcador radiopaco [PL] Znacznik radiocienicy [CS] Radioopákní značka [DA] Røntgenfast marker [FI] Radio-opaakki merkki [NO] Røntgentett marker [SV] Röntgentät markö [BG] [ET] Röntgenkontrastne märgistus [HU] Sugárfogó marker [LV] Rentgenkontrastains mar ieris [LT] Rentgeno spinduluose matoma žyma [RO] Marcaj radiopac [SL] Radioopabni označevalec
3. [EN] Balloon [ES] Balón [FR] Ballonnet [DE] Ballon [IT] Palloncino [TR] Balon [SK] Balón [PT] Balão [PL] Balon [CS] Balónek [DA] Ballon [FI] Pallo [NO] Ballong [SV] Ballon [BG] [ET] Balloon [HU] Ballon [LV] Balons [LT] Balion [SL] Balon
4. [EN] Stiffening wire [ES] Fiador [FR] Fil de rigidification [DE] Versteifungsdraht [IT] Filo di rinforzo [TR] Takviye teli [SK] Vysťuzovací drôt [PT] Fio de endurecimento [PL] Drut usztywniający [SV] Vytzužující drát [DA] Aftsvivende wire [FI] Jäikystyslanka [NO] Avatvingsvaier [SV] Förstärkningstrå [BG] [ET] Jäigastustraata [HU] Merevitézál [LV] Nosprieģošanas st ga [LT] Standžioji viela [RO] Fir de rigidizare [SL] Ojačitvena žica
5. [EN] Shaft [ES] Eje [FR] Gaine [DE] Schaft [IT] Asta [TR] Mil [SK] Telo [PT] Haste [PL] Trzon [CS] Dřk [DA] Skæft [FI] Varsi [NO] Skæft [SV] Skaf [BG] [ET] Vars [HU] Szár [LV] V rpstas [LT] Velenas [RO] Tijă [SL] Gred
6. [EN] Depth marking at 1015mm from the distal tip (TRP10015) [ES] Marcador de profundidad a 1015mm desde la punta distal (TRP10015) [FR] Marqueur de profondeur à 1015mm de l'extrémité distale (TRP10015) [DE] Tiefenmarkierung in einer Entfernung von 1015mm vom distalen Ende (TRP10015) [IT] Indicatore di profondità a 1015mm dalla punta distale (TRP10015) [TR] Distal uçtan 1015mm mesafedeki derinlik iüareti (TRP10015) [SK] Označenie hĺbky vo vzdialenosti 1015mm od distálneho konca (TRP10015) [PT] Marcação da profundidade a 1015mm da ponta distal (TRP10015) [PL] Oznaczenie głokości na 1015mm od

to narrow gaps and/or spaces and diminished safety and / or performance of the device and/or may lead to device failure, which in turn, may result in patient injury, illness or death.

- Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.
- Do not advance the TrapIT Trapping Balloon Catheter beyond the distal tip of the guide catheter.

Precautions

- The catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedure described do not represent ALL medically accepted protocols, nor are they intended as a substitute for the physician's experience and judgement in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.
- Handle with care to avoid product damage.
- Store in a dry, dark, cool place.
- All air must be removed from the balloon and displaced with contrast prior to inserting into the body.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon. Use only a contrast medium diluted to an appropriate concentration mixed solution of contrast medium and physiological saline. If the balloon cannot be positioned accurately there is a possibility that the OTW catheter will break if pulled out when the balloon is expanded over the distal end of the OTW type catheter.
- Before exchanging an OTW like catheter be sure all air is removed from the Hemostasis device/Y-connector
- Do not reinsert the catheter unless the balloon is fully deflated under vacuum. Apply negative pressure on the balloon to remove the liquid solution and remain negative pressure when inserting the balloon into the guide catheter. Moving the balloon in expanded state may damage the balloon or may break the shaft.
- Do not reinsert the TrapIT Trapping Balloon Catheter into the coil dispenser after procedural use.
- When using the 100 cm guide catheter compatible TrapIT in a 90 cm guide catheter use the shaft depth marker to position the TrapIT tip to avoid that the distal tip of the TrapIT is beyond the distal tip of the guide catheter entering the vasculature.
- When removing the protective tube from the balloon of this product, remove it carefully. Damage to the balloon may make it impossible or difficult for expansion and contraction of the balloon.
- Balloon pressure should not exceed the rated burst pressure (RBP 12 atm.).
- If resistance is felt, determine the cause before proceeding. Continuing to advance or retract the catheter while under resistance may result in damage / separation of the catheter.
- Exercise care during handling of the TrapIT Trapping Balloon Catheter to reduce the possibility of accidental breakage, bending or kinking of the device shaft.
- Do not twist or rotate the catheter more than one (1) full turn; the guide wire might be tangled.
- Do not use this product with two or more catheters in the guiding catheter at the same time.
- Use of a pressure-monitoring device is recommended to prevent over pressurization.

Reporting of adverse events

The user and / or patient should report any serious incident that has occurred in relation to this TrapIT to the manufacturer and the competent authority of state in which the user and / or patient is established.

Potential side effects / Adverse events

Possible side effects / adverse events include, but are not limited to the following

- Acute myocardial infarction
- Angina or unstable angina
- Arrhythmias, including ventricular fibrillation.
- Coronary vessel dissection, perforation, rupture or injury
- Drug reactions, allergic reaction to contrast medium.
- Embolism
- Infection
- Pyrogenic reaction
- Arterial Thrombosis
- Coronary artery spasm
- Death
- Hemorrhage or hematoma.
- Hypo / hypertension.

Package contains

One (1) TrapIT Trapping Balloon Catheter

Other items required but not provided

Interventional devices used for the intended purpose of the TrapIT are described in this IFU. The use and compatibility of the interventional devices with the TrapIT is described in the instructions supplied with the respective devices.

Preparation and inspection procedure

Note: Do not expose the catheter to organic solvents (e.g. alcohol).

Note: Store in a cool, dark and dry place

Warning: Do not use if the inner package is open or damaged.

1. Prior to use, carefully inspect the TrapIT catheter packaging and components for damage.
 2. Using sterile technique, transfer the dispenser coil with the TrapIT catheter into the sterile field.
- Complete the following steps to prepare the TrapIT for use
3. Gently detach the hub of the catheter from the packaging clip.
 4. Gently remove the catheter from the packaging hoop.
 5. Prepare an inflation device with the recommended contrast medium according to the manufacturer's instructions.
 6. Attach a stopcock to the hub balloon inflation lumen.
 7. Evacuate air from the balloon segment using the following procedure
 - a. Fill a 20 cc syringe or the inflation device with approximately 4 cc of the recommended contrast medium.
 - b. After attaching the syringe or inflation device to the stopcock on the hub, orient the balloon catheter with the distal tip and the balloon pointing in a downward vertical position.
 - c. Apply negative pressure and aspirate for 10 seconds. Slowly release the pressure to neutral, allowing contrast to fill the shaft of the balloon catheter.
 - d. Close the stopcock
 - e. Disconnect the syringe or inflation device from the inflation port of the balloon catheter.
 - f. Remove all air from the syringe or inflation device barrel.
 - g. Reconnect the syringe or inflation device to the stopcock on the inflation port of the balloon catheter.
 - h. Open the stopcock and maintain negative pressure on the balloon until air no longer comes into the syringe or inflation device.
 - i. Slowly release the device pressure to neutral.
 - j. Close the stopcock and disconnect the 20 cc syringe (if used) and connect the inflation device to the inflation port of the balloon catheter without introducing air into the system.

Warning: All air must be removed from the balloon and displaced with contrast prior to inserting into the body (repeat steps a through j, if necessary); otherwise, complications may occur.

Instructions for use

Note: Before insertion of the TrapIT into the guide catheter check size of the guide catheter and size of the OTW catheter to be removed or exchanged.

Guiding Catheter size	Maximum outer diameter of catheter to be removed or inserted
6F: Inner diameter of 0.070" or 1.78mm 3.3 French (1.10mm)	3.3 French (1.10mm)
7F: Inner diameter of 0.078" or 1.98mm 3.8 French (1.28mm)	3.8 French (1.28mm)
8F: Inner diameter of 0.088" or 2.24mm 4.6 French (1.56mm)	4.6 French (1.56mm)

1. While holding the guide wire in place, pull out the OTW catheter to be removed until the hub of the OTW catheter comes near the end of the guide wire.

Note: When inserting the TrapIT balloon catheter into the hemostasis device and guiding catheter, one's hand should support the balloon catheter and firmly grasp the proximal shaft. Shaft diameter differences should be taken into consideration when opening and tightening the hemostatic valve and upon withdrawal of the balloon catheter.

2. While holding the OTW catheter in position, insert the TrapIT through the Hemostasis device/Y-connector until the hard stop runs on the Hemostasis device/Y-connector.

Note: If a 100 cm compatible TrapIT is used in a 90 cm guide catheter, insert the TrapIT till the shaft depth marker is at the level of the Hemostasis/Y-connector cap.

3. Conform with fluoroscopy that the position of the proximal end of the TrapIT balloon is distally of the distal tip of the OTW catheter.

Note: It is important that the hemostatic valve is closed tightly enough to prevent blood leakage around the catheter shafts, yet not so tight that it restricts the flow of contrast into and out of the balloon or restricts catheter movement.

4. Inflate the TrapIT balloon to the recommended nominal pressure of 0.8 MPa (8 atm).
5. Pull lightly on the proximal end of the guide wire to conform the guide wire is fixed. If the guide wire is not secure enough, apply pressure until the wire is fixed. However, the pressure should not exceed the maximum expansion pressure of 1.2 MPa (12 atm)
6. Gently retract the OTW catheter while observing the position of the guide wire with fluoroscopy.
7. Tighten the hemostatic valve.
8. If it is still necessary to hold the guide wire for insertion of another OTW type catheter, back load the OTW type catheter over the guidewire and insert slowly to the vicinity of the balloon proximal part.
9. Deflate the TrapIT balloon till all liquid solution is removed from the balloon. Confirm deflation with fluoroscopy.
10. While maintaining the position of the guide wire carefully retract the TrapIT from the guide catheter and Hemostasis device/y-connector.
11. Tighten the hemostatic valve.

Note: After the deflated TrapIT balloon catheter is withdrawn, it should be wiped clean with gauze soaked with sterile, heparinized normal saline and stored. Prior to reinsertion, the balloon should be submerged in sterile, heparinized normal saline.

Note: If resistance is felt during removal of the balloon catheter from the guide catheter, the TrapIT balloon catheter, guide catheter, guide wire and OTW catheter should be removed together as a unit and set aside. Neither the TrapIT balloon catheter nor the guide catheter, guide wire or OTW catheter should be re-used for additional attempts, either separately or in combination.

12. Coil the balloon catheter after use and store.

Note: Care should be taken not to kink or bend the shaft upon coiling

Product disposal.

After use, this product may be a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

References
The physician should consult recent literature on current medical practice.

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[EN] Disclaimer of warranty and limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the IMDS product(s) described in this publication. Under no circumstances shall IMDS be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind IMDS to any representation or warranty except as specifically set forth herein.

Description or specifications in IMDS printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. IMDS will not be responsible for any direct, incidental or consequential damages resulting from the reuse of the product.

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IMDS Operations B.V.