

Polux™ SC 0.014" system

Excellent performance in routine
and challenging cases

Designed to treat resistant lesions

For a wide range of
indications

Semi-Compliant OTW Designed to treat BTK lesions

Optimized crossing profile

Best in class crossing capability for
diffused and challenging lesions

Slider™ hydrophilic coating


Ensures smooth entry to tight lesions

Kink resistant shaft

Smooth transmission through
tortuous vasculature

Low entry profile 0.017"

For reaching distal vasculature



Specifications

Balloon Diameter (mm)	Balloon Length (mm)												Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	5	10	15	20	30	40	60	80	100	120	150	200				
1.5	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
2.0	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
2.5	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
3.0	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
3.5	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
4.0	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
4.5	V	V	V	v	v	v	v	V	V	V	V	V	70	5	6	14
5.0	V	V	V	v	v	v	v	V	V	V	V	V	70	5	6	14
5.5	V	V	V	v	v	v	v	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14



BrosMed Medical Co., Ltd.

Add: 15th building, SMEs Park, Songshan Lake Hi-Tech Zone, 523808, Dongguan, Guangdong, China
 Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
 E-mail: sales@brosmed.com www.brosmed.com

BrosMed Medical B.V. (International Sales Office)

Mgr. Buckxstraat 8, 6134 AP Sittard The Netherlands
 Tel: +31 (0) 850 140 901
 E-mail: cs@brosmed.com
 Web: www.brosmed.com

US FDA K160256



GRA-6023/Rev.01